

BFH105

Inline Filter Holder



BFH105 inline filter holders are designed for fiber coupling, and ease of handling. It holds three standard $\text{Ø} 1 \times 5\text{mm}$ filters. The holder provides SMA 905 fiber coupler and is ideal for use with BTC100 series CCD spectrometer and BPS100 series constant photo intensity tungsten light sources.

Highlights

- Fiber Coupled
- Holds three $\text{Ø} 1 \times 5\text{mm}$ filters
- Can be used under ambient light

Specifications

Number of filters:	Three
Fiber coupler:	SMA905
Filter size:	$\text{Ø} 1 \times 5\text{mm}$
Optics wavelength range:	200 - 2400nm
Transmission:	> 60% without filter
Beam spot:	$\text{Ø} 5.5\text{mm}$

