



19 Shea Way, Newark, DE 19713
Tel: (302) 368-7824 Fax: (302) 368-7830
<http://www.bwtek.com> • E-mail: <http://www.bwtek.com/contact>

New Product – Econic: Scientific Stimulus Package

Company: **B&W Tek, Inc.**

Address: 19 Shea Way
Newark, DE 19713 USA

Telephone: +1-302-368-7824
Fax: +1-302-368-7830
Website: www.bwtek.com

Contact person: Jason Hammock
Email: jasonh@bwtek.com

New product Description (100 words):

B&W Tek presents a Scientific Stimulus Package for academic and government teaching labs. The Econic series comprises of high performance CCD spectrometers, each for less than \$900. Employing a 2000-element linear CCD array, a built-in 16-bit digitizer and a high speed USB 2.0 interface, this price can't be beat. Econic comes with a standard 90-day warranty (upgradeable). Standard wavelength ranges are 250-950 nm and 350-1050 nm with a standard slit size of 25 μm (upgradeable) delivering a resolution of 0.8-10 nm, depending on your configuration. External event synchronization and external trigger input are included, as well as BWSpec software. *econic.bwtek.com*