



19 Shea Way, Newark, DE 19713
Tel: (302) 368-7824 Fax: (302) 368-7830
<http://www.bwtek.com> • E-mail: <http://www.bwtek.com/contact/>

New Product – MiniRam II: Raman Spectrometer System

Company: **B&W Tek, Inc.**
Address: 19 Shea Way
Newark, DE 19713 USA

Telephone: +1-302-368-7824
Fax: +1-302-368-7830
Website: www.bwtek.com

Contact person: Jason Hammock
Email: jasonh@bwtek.com

New product Description (100 words):

B&W Tek, Inc is pleased to announce the next generation of portable Raman spectrometers - MiniRam II. This low cost, light weight battery operated patent pending Raman instrument leads the performance to cost ratio with a 10 cm^{-1} spectral resolution in $175 - 3100 \text{ cm}^{-1}$ Raman shift range, with a stabilized 785 nm excitation laser. By using a fiber coupled sampling probe and integrated handtop PC, a user can collect Raman spectra of solids and liquids at the scene where the information is most useful. Also available is an integrated cuvette based sampling system configuration. For more information, visit
