

B&W Tek, Inc information

Company: **B&W Tek, Inc.**

Address: 19 Shea Way
Newark, DE 19713

Telephone: 302-368-7824
Fax: 302-368-7830
Website: www.bwtek.com

Contact person: Jason Hammock
Email: jasonh@bwtek.com

Product Description:

BWID is specialized application software for all portable and desktop Raman spectrometers from B&W Tek. BWID and Raman spectroscopy are capable of recognizing materials at the molecular level. In practice, Raman spectra can be acquired through transparent packaging materials. BWID supports the identification of unknown samples as well as the verification of a known sample. The program quickly analyzes the substance in question and provides an immediate decision on its identification or validity. BWID features an intuitive user interface and managed workflow to minimize user error and ensure even novice users are quickly up and running. BWID assists user-definable methods for sample identification and user-built data libraries. The predefined method allows all instrument operators to perform sample identification through a one-click interface. Users can build unique spectral libraries to verify known materials or access third party libraries to identify unknowns. 21 CFR Part 11 compliance for pharmaceutical manufacturing.
www.bwtek.com/product/software/