

**Product Enhancement – Quest™ X CCD Spectrometer
Cylindrical Lens Option Now Available**

Company: **B&W Tek, Inc.**

Address: 19 Shea Way
Newark, DE 19713 USA

Telephone: +1-302-368-7824
Fax: +1-302-368-7830
Website: www.bwtek.com

Contact person: Jason Hammock
Email: jasonh@bwtek.com

New product Description (100 words):

Quest™ X is a low cost, high performance CCD spectrometer designed for UV-Vis-NIR applications with resolutions of 0.5 - 20 nm, depending on your configuration. Quest™ X employs a high speed USB 2.0 interface, a low stray light optical bench and an external trigger input. B&W Tek, Inc. is proud to announce Quest™ X has been enhanced to accommodate an optional cylindrical lens which provides a threefold increase in the spectral signal for a given integration time. With new dimensions of 102 x 67 x 34 mm, Quest™ X is bigger, better and brighter! For more information visit www.bwtek.com/product/spectrometer/questx.html

