

## New Product – ínnoram™: High Resolution Raman Spectrometer System

Company: **B&W Tek, Inc.**  
Address: 19 Shea Way  
Newark, DE 19713 USA

Telephone: +1-302-368-7824  
Fax: +1-302-368-7830  
Website: [www.bwtek.com](http://www.bwtek.com)

Contact person: Jason Hammock  
Email: [jasonh@bwtek.com](mailto:jasonh@bwtek.com)

### New product Description (80 words):

The ínnoram™ is the first portable Raman spectrometer system featuring an integrated touch-screen computer with research grade performance. It combines our core technology of high throughput optical design and deep TE-cooled CCD, delivering the performance of high sensitivity, extraordinary dynamic range, and spectral resolution of  $2\text{ cm}^{-1}$ . Together with our patented CleanLaze® laser and Raman probe with near-excitation cut-on at  $65\text{ cm}^{-1}$ , this system will equip you with the tools needed to tackle the dynamic Raman spectroscopy industry. Website: <http://www.bwtek.com/product/raman/innoram.html>

